




## Lost Connection Checklist

- 1 Check AC Adaptor is plugged into power outlet.
- 2 Power outlet is switched on.
- 3 USB Connecting Cable is connected to Wi-Fi Module (USB socket).
- 4 Power cable is connected to Wi-Fi Module (DC socket).
- 5 USB cable is connected to Bed Sensor Pad (BSP) Electronics Module.
- 6 Cycle the power - turn off the power to the Wi-Fi Module for at least 10 seconds and turn on again.
- 7 On the Wi-Fi Module, the white light should be permanently on and the green light should blink slowly.




### Further resources:

- Connection Diagram (overleaf).
- Lost Connection Checklist video is available here:  
[ensign.cura1.com/video/connection-checklist.mp4](https://ensign.cura1.com/video/connection-checklist.mp4)

**If the above does not resolve the lost connection, please contact your maintenance personnel.**

# Connection Diagram



## Lost Connection Checklist

- 1 Check AC Adaptor is plugged into power outlet.
- 2 Power outlet is switched on.
- 3 USB Connecting Cable is connected to Wi-Fi Module.
- 4 Power cable is connected to Wi-Fi Module (DC socket).
- 5 USB Connecting Cable is connected to Bed Sensor Pad (BSP).
- 6 Cycle the power - turn off the power to the Wi-Fi Module for at least 10 seconds and turn on again.
- 7 On the Wi-Fi Module, the white light should be permanently on and the green light should blink slowly.